

F

WOUND DRESSING PRODUCTS

F
39

WOUND DRESSING PRODUCTS

F001 Gauze Swabs

Material: 100% Cotton
Cotton Yarn: 40'sx40's
Mesh: 18x11, 18x14, 19x9, 19x11, 19x15, 20x12, 24x20, 26x18, 28x24 or other mesh per customer's demand.
Size: 2"x2", 3"x3", 4"x4", 4"x8" or other size per customer's demand.
Ply: 4Ply, 8Ply, 12Ply, 16Ply, etc.
With or Without X-Ray
Edge: Folded or unfolded
Color: White, Blue, Green
Sterile or Non-Sterile
Sterilization: EO or Gamma
Sterile Packing: 1pc/Two sides of paper pouch or one side paper+one side transparent PE film pouch, 2pcs/Pouch, 5pcs/Pouch or other packing per customer's demand
Non-Sterile Packing: 100pcs/Paper pack or Polybag



F002 Gauze Roll

Material: 100% Cotton
Cotton Yarn: 40'sx40's
Color: White Bleached, Chlorine Free, Non-Toxic, Non Abrasive, Latex Free, Clean And Dust Free
Mesh: 18x11, 18x14, 19x9, 19x11, 19x15, 20x12, 24x20, 26x18, 28x24 or other mesh per customer's demand.
Size: 36"x50m, 36"x100yards, 36"x100m, etc.
Ply: 1Ply, 2Ply, 4Ply etc.
With or Without X-Ray
Packing: 1roll/blue kraft paper pack, then each roll in nylon bag



F003 Laparotomy Sponge

Material: 100% Cotton
Cotton Yarn: 40'sx40's
Color: White Bleached, Chlorine Free, Non-Toxic, Non Abrasive, Latex Free, Clean And Dust Free
With X-ray Detectable or Without X-ray Detectable
Mesh: 18x11, 18x14, 19x9, 19x11, 19x15, 20x12, 24x20, 26x18, 28x24 or other mesh per customer's demand.
Size: 22.5x22.5cm, 30x30cm, 45x45cm, etc, or other size customized
Ply: 4Ply, 8Ply, 12Ply etc.
Sterilization: EO or Gamma, Sterile or Non-Sterile
Sterile Packing: 5Pcs/Two sides of paper pouch or one side paper+one side transparent PE film pouch.
Non-Sterile Packing: 20pcs/Polybag or in Bulk




F045 Paraffin Impregnated Gauze Pad

Description features: Sterile Paraffin Gauze Dressings.
 Are made of 100% cotton gauze swabs impregnated with white soft paraffin.
 Their multi-trulle structures are indicated as a primary wound contact dressing for manufactured from high-absorbent gauze swabs which have folded edges and thus prevent loose fibres and shredding.
Size: 10cmx10cm etc.
Packing: Each piece in foil pouch, 10pouch/Inner box, 100pcs/ctn
Carton Size: 62x28x32cm G.W: 8.0kgs N.W: 7.0kgs



F005 Wow Gauze Bandage

Material: 100% Cotton
Cotton Yarn: 40'sx40's
Color: White Bleached, Chlorine Free, Non-Toxic, Non Abrasive, Latex Free, Clean And Dust Free
Mesh: 18x11, 18x14, 19x9, 19x11, 19x15, 20x12, 24x20, 26x18, 28x24, or other mesh per customer's demand
Width Size: 1", 2", 3", 4", 5", 6"
Length Size: 3yds, 4yds, 5yds, 6yds, 10yds or other size
Sterilization: EO or Gamma, Sterile or Non-Sterile
Sterile Packing: Each roll in transparent pe airproof pouch or other packing per customer's demand
Non-Sterile Packing: 1doz/Blue kraft paper packing



F006 Non-Woven Swabs

Introduction: It is made from 30% terylene and 70% glue all cut edges are folded or unfolded and stitched evenly with a steady loop attached firmly to one corner. They are most useful for bleeding control or tissue protection in the operation rooms. It is normally with/without X-ray detectable threads or tapes and can be dyed in various colors. It meets USP and European Pharmacopoeia

Gram of Material: 30g/m², 40g/m², etc.

Size: 2"x2", 3"x3", 4"x4", 4"x8" or Other Size Per Customer's Demand.

Ply: 4Ply, 8Ply, etc.

Color: White, Blue, Green

Sterile or Non-Sterile

Sterilization: EO or Gamma

Sterile Packing: 1pc/Two sides of paper pouch or one side paper+one side transparent pe film pouch, 2pcs/pouch, 5pcs/pouch or other packing per customer's demand.

Non-Sterile Packing: 100pcs/Paper Pack or 100 pcs/Polybag



F007 Absorbent Cotton Wool

Introduction: Made from 100% natural cotton after cleaning purifying and bleaching, which is used for cleaning the wound. It can be made into roll.

Grade of Cotton: Natural or super absorbent

Weight/Roll: 50g, 100g, 125g, 150g, 200g, 400g, 454g, 500g, 1000g, etc.

Packing: Each roll in Blue kraft paper packing

Specification	Roll/Blue kraft paper	QTY	Carton Size (cm)
50g	1	500	58x52x51
100g	1	200	61x31.5x61.5
200g	1	100	62x43.5x40
250g	1	100	61x49x49
400g	1	36	60x38x46
454g	1	36	57x31x57
500g	1	36	49x31x49
1000g	1	25	63x31x63



F010 Gauze Ball

Specification: 40's 18x14mesh or 19x11mesh or 19x15mesh or 24x20mesh or 26x18mesh etc.

Size: Dia.3cm or Dia.4cm or Dia.5cm are available

Shape: Ball shape or peanut shape are available

Color: Pure white, green, blue are available



F011 Cotton Swab

Material Of Stick: PP, PS, wooden

Type: Single head, double head

Diameter of The Plastic Stick: 2.5mm-4.0mm

Length of The Plastic Stick: 75mm-400mm

Diameter of The Cotton Tip: 4-15mm

Length of The Cotton Tip: 15-50mm

3"	100pcs/bag	10bags/box	20boxes/ctn	47x37x23cm
6"	100pcs/bag	10bags/box	20boxes/ctn	43x37x37cm
6"	50pcs/bag	20bags/box	8boxes/ctn	63x37x49cm
8"	50pcs/bag	20bags/box	8boxes/ctn	63x37x58cm



F038 Self-Adhesive Eye Pad

Material: Non-woven

Color: White, Flesh

Weight Of Adhesive: 30-32g/m²

Shape: Oval

Size: 80x50mm, 85x55mm, 95x65mm, etc.

Packing: 50pcs/box, 24box/ctn

Carton Size: 51.5x27.5x24cm



F008 Zig Zagcotton Wool

Introduction: Made from 100% natural cotton after cleaning purifying and bleaching, which is used for cleaning the wound. It can be made into piece.

Grade of Cotton: Natural or super absorbent

Weight/Roll: 25g, 50g, 100g, 200g, 250g, 500g, 1000g, etc.

Packing: Sealed poly bag or drawstring bag, self-sealing bag

Specification	Bags/ctn	Carton Size (cm)
25g	400	58x52x51
50g	200	64x58x46
100g	100	63x46x46
200g	50	63x46x48
250g	30	70x42x48
500g	20	46x40x47
1000g	10	68x50x43



F009 Cotton Ball

Introduction: Made from 100% natural cotton after cleaning purifying and bleaching, which is used for cleaning the wound.

Grade of Cotton: Natural or Super Absorbent

Weight/Roll: 0.3g, 0.5g, 1.0g, 2.0g, 3.0g, 4.0g, 5.0g, etc.

Sterilization: EO or Gamma, Sterile or Non-Sterile

Packing: Sealed poly bag packing

Specification	Piece/bag	Bags/ctn	Carton Size (cm)
0.3g/pc (Sterile)	50	660	64x58x48
0.3g/pc (Sterile)	100	330	64x58x46
0.3g/pc (Non-sterile)	1000	33	64x58x46
0.3g/pc (Non-sterile)	2000	16	64x58x46
0.5g/pc (Sterile)	50	400	64x58x46
0.5g/pc (Sterile)	100	200	64x58x46
0.5g/pc (Non-sterile)	1000	20	64x58x46
0.5g/pc Non-sterile)	2000	10	64x58x46
2g/pcs (Sterile)	50	100	64x58x46
2g/pcs (Sterile)	100	50	64x58x46
2g/pc (Non-sterile)	1000	5	64x58x46
2g/pc (Non-sterile)	2000	1	46x40x38



F015 Spandex Elastic Bandage

Material: 85% cotton With 15% spandex Crepe or Plain

Color: natural, Flesh or Bleach

Weight: 75g/m², 80g/m², 85g/m², 90g/m², etc.

Width: 2", 3", 4", 5", 6" etc.

Length: 4m, 4.5m, 5m, 6m, etc.

Sterile or Non-Sterile

Packing: 1roll/Polybag, 12roll/bag



F016 PBT Elastic Bandage

Material: 65% cotton with 35% polyester

Color: Bleach

Weight: 28g/m², 29g/m², 30g/m² etc.

Width: 2", 3", 4", 5", 6" etc.

Length: 4m, 4.5m, 5m, 6m, etc.

Sterile or Non-Sterile

Packing: 1roll/Polybag, 12roll/bag



F017 Cotton Elastic Bandage

Material: 100% cotton
Color: Bleach
Weight: 75g/m², 80g/m², 85g/m², 90g/m², etc.
Width: 2", 3", 4", 5", 6" etc.
Length: 4m, 4.5m, 5m, 6m, etc.
Sterile or Non-Sterile
Packing: 1roll/Polybag, 12roll/bag



F018 High Elastic Bandage

Material: 62% cotton With 38% polyester
Color: Bleach, Natural
Weight: 80g/m², 100g/m², etc.
Width: 2", 3", 4", 5", 6" etc. Length: 4m, 4.5m, 5m, 6m, etc.
Sterile or Non-Sterile
Packing: 1roll/Polybag, 12roll/bag



F071 Fixing Tape Roll

Material: Non-Woven or PU

Size	roll/box	boxes/ctn	Carton Size
5cmx10m	1	160	61.5x39x23cm
10cmx10m	1	80	39x31.5x43cm
15cmx10m	1	60	39x31.5x47.5cm
20cmx10m	1	40	39x31.5x42.5cm



F046 First-Aid Bandage

Cotton Filled Size: L, M, S
Packing: Vacuumize
Packing:
Small 1pc/Vacuumize pack, 400pack/ctn 50x44x36cm
Large 1pc/Vacuumize pack, 400pack/ctn 50x46x43.5cm



F019 Net Bandage

Material: Acrylic with elastic
Color: White
Weight: 120g/m², etc.
Size: 0.8cmx25m, 1.7cmx25m, 2cmx25m, 2.2cmx25m, 2.6cmx25m, 3.8cmx25m, 4.4cmx25m, 5.5cmx25m, 6.8cmx25m, 8cmx25m



F020 Non-Woven Self-Adhesive Bandage

Material: Non-woven PP
Color: White, blue green, red, brown, yellow etc.
Width: 2", 3", 4", 5", 6" etc. Length: 4m, 4.5m, 5m, 6m, etc.
Sterile or Non-Sterile
Packing: 1roll/Polybag, 12roll/box



F033 Ice Bandage

1. Ice wrap is a cold bandage which uses a revolutionarily new technique for treating injuries ,muscle pain and strains.
2. This bandage reduces recovery time from injuries by combing both the cooling and compressing process in an easy-to-use way.
3. A specially formulated solution is infused into the elastic bandage, then begins to evaporate upon contacting skin, thereby causes cooling of the applied area.
4. Icewrap reduces swelling and provides sufficient cooling to make it a very effective treatment.

Packing:

10cmx3.2m	120rolls/ctn	66x33x26cm
10cmx3.2m	48rolls/ctn	36x27x42cm
7.5cmx3.2m	120rolls/ctn	66x33x26cm
7.5cmx3.2m	80rolls/ctn	36x27x42cm



F024 Trauma Wound Dressing

Feature: High hook strength, good elasticity, high water absorption
Material: Cotton, Ployester, Plastic
Color: Green/OEM
Size: 10cm/15cm/20cmx4.5m
Weight: 90-120g
Package: Double layer vaccum package,100pcs/ctn



F021 Cotton Self-Adhesive Bandage

Material: Cotton
Color: White,blue green, red, brown, yellow etc.
Weight: 75g/m², 80g/m², 85g/m², etc.
Width: 2", 3", 4", 5", 6", etc.
Length: 4m, 4.5m, 5m, 6m, Etc.
Sterile or Non-Sterile
Packing: 1roll/Polybag, 12roll/box



F022 Plaster Of Paris Bandage

Material: Cotton with plaster
Color: White
Weight: 390g/m², 400g/m², 410g/m², etc.
Width: 2", 3", 4", 5", 6", etc. Length: 4m, 4.5m, 5m, 6m, etc.
Sterile or Non-Sterile
Packing: 1roll/poly Bag, 12roll/box



F035 Fiberglass Casting Tape

Material: Fiberglass
Size: 5x360cm, 7.5x360cm, 10x360cm, 12.5x360cm, 15x360cm
Packing: 1pc/Aluminum foil bag, 50pcs/box, 100boxes/ctn




F032 Sports Tape

Ideal for sport protection, water proof, permeability
Size: 2.5cmx10y, 3.8cmx10y, 5cmx10y, 2.5cmx15y, 3.8mx15y, 5cmx15y
Packing: 6rolls/box, 144boxes/38x30.5x34cm
1roll/box, 125rolls/38x38x32.3cm



F074 Orthopedic Wool Padding

Material: 100% Cotton or synthetic fiber
Size: 5cm, 7.5cm, 10cm,15cm
Length: 2.7m, 3.0m, 4.0m



F048 Triangle Bandage

Material: 100% Cotton or rayon fabric
Cotton Yarn: 40x34, 50x30, 48x48, etc.
Size: 36"x36"x51", 40"x40"x56"
With or without safety pin colored: Grey or Bleached.
Packing: 1pc/Sealed poly bag, 240bags/ctn
Carton Size: 45x38x22cm, 6.5/5.0kgs



F025 Zinc Oxide Plaster

Material: 10% Zinc oxide with 90% hot melt adhesive
Color: White, Flesh.
Weight of adhesive: 80-85g/m²
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm, etc.
Length: 5m, 10m, etc.
Packing: Plastic can, tinplate or simple packing



F026 Non-Woven Surgical Tape

Material: Non-woven with adhesive
Color: White, Flesh
Weight of adhesive: 30-32g/m²
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm, etc.
Length: 5m, 10m, etc.
Packing: Plastic can, plastic cutter or simple packing



F027 PE Surgical Tape

Material: Polyethylene with adhesive
Color: Transparent
Weight of adhesive: 30-34g/m²
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm, etc.
Length: 5m, 5y, 10y, etc.
Packing: Plastic can, Plastic cutter or simple packing.




F028 Silk Surgical Tape

Material: Silk with adhesive
 Drilled or non-drilled is available
 Line shape or zig shape is available
Color: White
Weight of adhesive: drilled 50-54g/m², Non-drilled 32-34g/m²
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm, etc.
Length: 5m, 10m, etc.
Packing: Plastic can, simple packing.




F029 Wound Adhesive Plaster

Material: PVC, PE, PU, Cotton, Elastic, Non-woven etc.
Color: White, Flesh, etc.
Weight of adhesive: 32-40g/m²
Shape: Round, Square, Rectangle, Butterfly etc
Size: 22x22mm, 25x25mm, 40x10mm, 38x38mm, 60x19mm, 72x19mm, 76x19mm, 76x38mm, 114x50mm, etc.
Packing: 50pcs/Box or 100pcs/Box, 100boxes/ctn



F030 Drilled Plaster

Material: Drilled Zinc oxide with hot melt adhesive.
 Waterproof or non-water proof, With drug or without.
Color: White, Flesh
Weight of adhesive: 60-65g/m²
Width: 2.5cm, 5cm, 7.5cm, 10cm, 12cm, etc.
Length: 1m, 5m, etc.
Packing: 1roll/Box, 80roll/ctn



F031 Healing Capsicum Plaster

Capsicum plaster uses the full cotton cloth as the raw material, which was coated on capsicum consisted zinc oxide adhesive. Capsicum plaster can make the partial skin generate physical reaction to accelerate blood fluxion
 This plaster has drilled type and not-drilled type.
 All the products adopt the anti-sticky paper and film as the interleaver.
 Can be used by one piece or cut size conveniently and easily.




F037 Infusion Plaster

Material: Non-woven
Color: White, Flesh
Weight Of Adhesive: 30-32g/m²
Shape: Rectangle
Size: 75x25mm, 75x35mm, etc.
Packing: 5pcs/pack, 20packs/box, 100boxes/ctn
Carton Size: 42x39.5x37cm



F1038 IV Cannula Dressing

Material: Transparent film dressing
 Hydrocolloid I.V dressing
 PE Foam I.V dressing
 Non-woven I.V dressing
Packing: 1pc/bag, 50pcs/box, 40boxes/ctn



F036 Self-Adhesive Wound Dressing

Material: PU, Non-woven, etc.
Color: White, Flesh, etc.
Weight of adhesive: 70-75g/m²
Shape: Round corner, Square corner.
 With mat or without is available.
Size: 5cmx7cm, 6cmx8cm, 10cmx10cm, 10cmx12cm, 10cmx15cm, 10cmx20cm, 10cmx25cm, 10cmx30cm, 10cmx35cm, etc.
Packing: 50pcs/box



F013 Cosmetic Cotton Swabs

Introduction: Made from 100% Super Absorbent cotton and 100% natural paper fibre, the material is very soft. Purifying with high temperature, which is used for cleaning the wound or dressing.
Shape: Round or square
Type: Plain, lattice and flower pattern Concavo-convex dot, etc.
Size: 5.5cmx5.5cm, 6cmx6cm, 7.5cmx5cm
Packing: 80pcs/Sealed plastic bag, 100pcs/Sealed plastic bag, 100pcs/box or other package per customer's demand
Normal Packing:
 Size 6.0cm(Round) 100pcs/Plastic bag 100bags/ctn 62x32x38cm
 Size 6.0x6.0cm(Square) 100pcs/Plastic bag 100bags/ctn 62x32x38cm



F039 Alcohol Prep Pad

Contents: 70% Isopropyl Alcohol, Non-woven
Size: 60x30mm, 65x30mm, 65x58mm, 140x180mm or other size customized
Type: 1Ply, 2Ply, 4Ply, etc.
Packing Size: 50x50mm, 60x85mm, etc.
Packing Quality: 1Pc/Aluminum foil bag, 200pcs/box, 20000pcs/ctn 53x45.5x29cm
 100pcs/box, 10000pcs/ctn 42.5x28.5x28.5cm



F040 Povidone-Iodine Prep Pad

Contents: 10% Povidone-iodine Solution Non-woven
Size: 65x30mm, 65x58mm, or other size customized
Type: 1Ply, 2Ply, 4Ply, etc.
Packing Size: 50x50mm, etc.
Packing: 1pc/Airtight Pack , 100pcs/box, 10000pcs/ctn
Carton Size: 54x24x29cm



F041 Adhesive Tape Remover Pad

Contents: Adhesive Tape Remover Pad
Non-Woven Size: 65x30mm
Packing Size: 50x50mm
Packing: 1pc/Airtight Pack , 100pcs/box, 10000pcs/ctn
Carton Size: 54x24x29cm



F042 Sting Relief

Sting relief used in wiping the skin area where was stung by mosquito, bee and other insects for relief of pain and itch.
30x60mm x2 non-woven pad wrapped in a foil paper sachet sized 50x50mm
Carton Size: 54x23x28cm



F043 Povidone-Iodine Swabstick

10% Povidone-iodine
Plastic Stick
Packing:
1pc/pack, 50pcs/box, 1000pcs/44x31x35cm
3pcs/pack, 75pcs/box, 1500pcs/44x31x35cm



F1039 Sponge Clean Stick

1# 18cm Flat Stick 20pcs/bag 40bags/ctn 53x52x41cm
2# 15cm Flat Stick 20pcs/bag, 40bags/ctn 44x34x42cm
3# 10cm Circl Stick 50pcs/bag, 20bags/ctn 47x16.5x40cm



F057 Basic Dressing Set1 (Paper-Plastics Compound Bag)

Packing: 100bags/ctn 47x31x41cm

1pc	Dressing Tray
5pcs	Cotton Wool Balls 0.5gm
5pcs	Gauze Swab 7.5cmx7.5cm-5pcs
3pcs	Dressing Forceps
1pc	Surgical Drape
1pc	Hand Towel
1pc	Limpet Bag
1pc	Sterile Field



F059 Basic Dressing Set3(Blister Bag)

1pc	Dressing Tray
6pcs	Gauze Ball
2pcs	Dressing Forceps
1pc	Surgical Drape
1pc	Hand Towel



F061 Dressing Set

Packing: 40bags/ctn 38x29x28cm

1pc	PE+PAPER 33x45cm
2pcs	Forceps 4#(blue and green)
5pcs	40g/m ² (D3cm)0.8G Non-woven ball



F058 Basic Dressing Set2(Shrink Bag)

Packing: 200bags/ctn 63x29x34cm

1pc	Dressing Tray
8pcs	Cotton Wool Balls 0.5gm
5pcs	Gauze Swab 7.5cmx7.5cm
3pcs	Dressing Forceps
1pc	Surgical Drape
1pc	Sterile Field



F060 Dressing Set

Packing: 120bags/ctn 41x34x34cm

1pc	Forceps 4#
1pc	Surgical blade
3pcs	Gauze swab 5x5cm



F062 Dressing Set

Packing: 40bags/ctn 38x29x28cm

1pc	PE+PAPER 33x45cm
2pcs	Forceps 4#(blue and green)
5pcs	20x20cm(28x24mesh)1.2G Gauze ball
5pcs	Non-woven swab 7.5cmx7.5cm



F063 Dressing Set

Packing: 200bags/ctn 69x42x34cm

2pcs	Forceps 4#(white and green)
3pcs	Gauze ball 28x24mesh(15x20cm)



F065 Dressing Set

Packing: 40bags/ctn 39x30x30cm

1pc	PE+PAPER 33x45cm
1pc	Hemostasia Scissor
1pc	Forceps 4
5pcs	20x20cm(28x24mesh)1.2G Gauze ball
5pcs	Non-woven swab 7.5cmx7.5cm



F069 First Aid Kit

Material: 600D oxford
Size: 32x22x7cm
Can be configured according to customer requirements
Packing: 20pcs/ctn 81x34x47cm

Item	Specification
ICE bag	160g
PVC Gloves	M, powdered, nature
Non-woven Surgical Tape	1.25cmx5y, white, simple packing
Non-woven Surgical Tape	5cmx5y, white, simple packing
Wound Adhesive Plaster	72x19mm, PVC, brown color
PBT Elastic Bandage	3" x4m, bleach, cellophane packing
PBT Elastic Bandage	5cmx4m, cellophane packing
Gauze Bandage	5cmx5m
Gauze Bandage	7.5cmx5m
Triangular Bandage	Mesh 46x30, 36"x36"x51"
Breathing Mask With One-Way Valve	20cmx20cm
Sterile Gauze Swab	Mesh 24x20, 10cmx10cm-12p, 5pcs/bag
Sterile Gauze Swab	Mesh 24x20, 7.5cmx7.5cm-12p,5pcs/bag
Applicator	3", 6pcs/bag
Cotton Wool	3pcs/roll
Alcohol Prep Pad	65x30mm-2ply
Scissor	/
Tweezer	4#



F066 Dressing Set

Packing: 40bags/ctn 39x30x30cm

1pc	PE+PAPER 33x45cm
1pc	Forceps 4
5pcs	20x20cm(28x24mesh)1.2G Gauze ball
5pcs	40g/m²(D3cm)0.8G Non-woven ball



F051 Surgical Hand Brush

Size: 8x5x4cm
Packing: 1pc/Blister bag, 250pcs/ctn 48x38x32cm



F070 First Aid Kit

Material: 420D oxford
Size: 21x13x6cm
Packing: 180pcs/ctn 62x46x30cm

Item	Specification
ICE bag	160g
Non-woven Surgical Tape	1.5cmx5y, white,simple packing
Wound Adhesive Plaster	72x19mm, PVC,brown color
PBT Elastic Bandage	3" x4m, bleach, cellophane packing
Gauze Bandage	5cmx5m
Gauze Bandage	7.5cmx5m
Triangular Bandage	Mesh 46x30, 36"x36"x51"
Breathing Mask With One-Way Valve	20cmx20cm
Sterile Gauze Swab	Mesh 24x20, 7.5cmx7.5cm-12p, 5pcs/bag
Applicator	3", 6pcs/bag
Cotton Wool	3pcs/roll
Alcohol Prep Pad	65x30mmx2ply
Scissor	/
Tweezer	4#



F052 Iodine Hand Brush

Size: 8x5x4cm
Specification: 7.5% iodine, with nail swabs
Package: 1pc/blister bag, 30pcs/box, 300pcs/ctn
Carton Size: 58x37.5x40cm



F053 Chlorhexidine Hand Brush

Size: 8x5x4cm
Specification: 4% chlorhexidine, with nail swabs
Package: 1pc/blister bag, 30pcs/box, 300pcs/ctn
Carton Size: 58x37.5x40cm



F068 Hemostasis Pincers

Size: Length 14cmxwidth7.5cm
Material: ABS and nylon are available
ABS for EO sterilization
Nylon material for high temperature steam sterilization.
Color blue, green, black, red, etc. are available
Packing: individual sterile Blister pack
MOQ: 20000pcs/color

F067 Plastic Forceps

Size	pcs/bag	bags/ctn	Carton Size
1# 10.8cm	200	25	64x46x46cm
2# 14cm	200	20	64x46x46cm
3# 11.5cm	200	10	51x32x35cm
4# 12.5cm	200	10	50x42x30cm
5# 18cm	1000	/	50x38x42cm
6# 10.8cm	200	25	64x46x46cm

