




G HOSPITAL UNIFORM

G001 Round Nurse Cap

Material: Non-woven Polypropylene
Weight: 10g, 12g, 14g, 16g, 18g, etc.
Size: 18", 21" or 24"
Color: White, Green or Blue, Pink, Yellow, Red, etc
Sterile Method: None or EO Sterile
Packing: 100pcs/Poly bag, 10bags/ctn
Carton Size: 24x24x34cm



G002 Clip Cap

Material: Non-woven Polypropylene
Weight: 10g, 12g, 14g, 16g, etc.
Color: White, Green Blue, Pink, Red, Yellow
Size: 20", 21" or 24"
Type: Single elastic/double elastic
Sterile Method: None or EO Sterile
Packing: 100pcs/Polybag, 10bags/ctn
Carton size: 53x24x35cm



G004 Doctor Cap

Material: Non-woven Polypropylene
Nylon Net Fabric
Color: Blue
Type: Tie or elastic back
Size: 18", 21" or 24"
Packing: 100pcs/Polybag, 10bags/ctn
Carton Size: 51x29x53cm



G006 Nurse Paper Face Mask

Material: High quality wet strength paper
Comfortable elastic ear loops for easy donning product
Ply: 1ply, 2ply
Size: 20x7cm
Color: White
Packing: 100pcs/Box, 100boxes/ctn
Carton Size: 52x38x24cm



G003 Doctor Cap

Material: Non-woven Polypropylene
Weight: 20g, 25g, 30g, 40g, etc
Color: Green/Blue
Type: Tie or elastic back
Size: 61x14cm; 64x23cm etc
Latex free
The SMS, Spun Lace Fabric is available
Packing: 100pcs/Polybag, 10bags/ctn
Carton Size: 50x50x30cm



G005 Non Woven Face Mask

Material: Non-woven polypropylene
Conforms to american nelson n95 test standards.
Filtration efficiency of dust and powder is over 95%.
Type: Ear-Loop and With Tie
Ply: 1ply, 2ply or 3ply
Detailed Weight: 1Ply-38g/m², 2Ply-25g/m²+18g/m²
3Ply-25G/m², +20g/m²+18g/m²
Size: 18x9cm
Color: White, Green, Blue, Pink etc
Packing: 50pcs/Box, 40boxes/ctn
Carton Size: 42x39x37cm



G007 Dust-Proof Face Mask

Material: Wood and any toxic free protection non-woven pp
Conforms to FFP1 and FFP2. To the skin adjustable contoured aluminum.
Nosepiece helps to ensure proper fit.
Weight: 16g, 18g, 20g, etc.
Color: White, Blue, Dark Blue
Packing: 50pcs/Inner Box, 20boxes/ctn
Carton Size: 60x26x25cm



G1008 Surgical Hood

Material: Non-Woven PP
A disposable surgical hood which conforms to the wearer's head. Is usable either with a traditional surgical mask alone or with such a mask.
In combination with a laminar air flow face plate contains bacterial fallout from the wearer's head.
It helps maintain sterile conditions in hospital operating rooms.
Weight: 16g, 20g, 25g, 30g, etc.
Size: 30x24cm
Color: White, Green, Blue
Packing: 100pcs/Poly bag, 10bags/ctn
Carton Size: 38x30x30cm



G008 Face Guard

Material: Non-Woven PP
Conforms to american nelson N95 test standards.
Ear-Loop Type: 2ply or 3ply
Cloth, low breathing resistance and high filtration Efficiency. And PVC transparent shield.
Mini-fiber high electrostatic filter, filtration efficiency of dust and powder is over 95%.
Repel continually contamination of dangerous blood And saliva continually. Inner cover is well sweet-Absorbed, while outer cover is water resistant.
Packing: 25pcs/Inner box, 20boxes/ctn
Carton Size: 40x34x21cm



G010 N95 Face Mask(With/Without Valve)

Material: Polypropylene
Conforms to international standards
EN149: 2001 FFP1 FFP2 CE0194
AS/NZS 1716: 1994 NIOSH N95
Low breathing resistance with high filtration efficiency.
Adjustable contoured aluminum nosepiece and a and increase Worker comfort.
Easy-breathing exhalation valve with inner silicon-piece.
Can facilitate the air exchange and increase worker comfort.
Special 3-ply design respirator produced from high filtrable. And electrostatic material for resisting normal toxic and tiny micro particle.
Suggested Applications:
Spray-paint/printing, dyeing, garage, factories of electronic aid and food disinfecting and the insalubrious workshops.
Packing: 20pcs/Box, 20boxes/ctn
Carton Size: 62x26x54cm



G009 New Type Dust Proof Face Mask

Description:
Low breathing resistance with high filtration efficiency.
Adjustable nose seal, easy to fold and carry.
Special 3-ply design respirator produced from high filtrable and electrostatic material with high efficient activated carbon cloth for relief from the irritating effects of nuisance odor, acid gas, hydrogen fluoride and ozone.
Suggested Applications:
Polluted environment, especially good for absorbing and filtrating the exhausted gas from motor vehicles.
Packing: 500pcs/Inner box, 20boxes/ctn
Carton Size: 60x27x25cm

G011 Filtering Half Mask

Flat-folded CE EN 149: 2001
Protection against solid and fluid particles (FFP1 and FFP2)
And protection against toxic particles(FFP3 only), single packed without valve, FFP1/with valve FFP2/with valve, FFP3
Packing: 20pcs/Box, 400pcs/ctn
Carton Size: 57x43x32cm



G012 Surgeon Gown

Material: Non-Woven PP/PP+PE+PE/SMS
Features:
1. With ties at back of the neck and waist.
2. Comfortable sleeve tailoring for maximum freedom of movement.
3. Soft knit cuffs for wrist comfort, knitted cuff is available.
4. Secure sleeve seams for protection.
5. Extra protection reinforcement liner on sleeve and chest and full-body is available.
6. Vecro neck back closure, side waist band.
7. 1 hand towel is available.
8. Anti-bacterial alcohol-repellent, blood-repellent waterproof dustproof, breathable.
Color: White/blue/Yellow/Green
Weight: 40/45/50GSM PP or SMS
Size: S 110x154cm, M 120x150cm, L 130x155cm, XL 140x160cm
Sterile Method: Non Sterile or EO Sterile
Packing: 1pc/Blister bag, 50pcs/ctn



G013 Uniform Isolation Gown

Material: Non-Woven PP/PP+PE
Full length for maximum coverage, Latex Free.
Sewn seams provide superior strength.
And extra length on waist ties allow gown to be easily secured in front.
With ties at back of the neck and waist.
Elastic wrist or with knitted cuff.
Suited to lab, industry, hospital and daily home use.
Color: White, Blue, Yellow, Green
Weight: 18/25/30/40GSM or 18/25/30/40 GSM PP+15PE
Size: 115x137cm or 120x140cm etc
Sterile Method: EO Sterile
Packing: 10pcs/Poly bag, 100pcs/ctn or 1pc/ Sterile packing, 50pcs/ctn
Carton Size: 56x32x35cm



G015 CPE Gown

Material: CPE
Description:
Prevent and isolate dust, particles, Alcohol, blood, bacteria and virus invading. Can be used as a visit coat. Prevent the cross infection in medical treatment at sanitation field. Accord with the hygienic standard. Characteristics: prevent and isolate dust, particle, alcohol, blood bacteria and virus invading. Thumb loop wrists prevent sleeves from sliding up when putting gloves on.
Color: Blue
Size: 125x200cm, 62g
Packing: 1pc/bag, 100pcs/ctn
Carton Size: 46x20x25cm



G016 Lab Coat

Material: Non-Woven PP/PE
Color: White, blue, Yellow, Green
Weight: 30/35/40/45GSM
Size: M 105x130cm, L 105x135cm, XL 115x140cm, XXL 115x145cm
 With label collar/knitted collar and 4 press buttons button.
 Closure or zipper closure or velcro closure.
 Suited to hospital and daily home use.
Packing: 10pcs/Poly bag, 100boxes/ctn
Carton Size: 56x32x42cm



G018 Disposable Scrub Coat

Material: Non-Woven PP or SMS material
 Disposable scrub coat are attractive as well
 Suited to hospital and daily home use.
 V collar top
 Elastic waist pant
Fabric: 48g, 55g
Color: Green, blue
Packing: 1pcs/Inner polybag, 50boxes/ctn
Carton Size: 43x34x55cm



G071 Sterile Delivery Pack

1pc Non-woven surgical gown
 1pc 86x80cm Packing clothing 2ply(Non-woven, thin film)
 1pc 97x73cm Sheet 3ply(thin film, Non-woven, paper pulp)
 2pcs 97x73cm Triangle trouser 3ply(thin film, Non-woven, paper pulp)
 1pc 66x43cm Non-woven cloth 2ply(thin film and Non-woven)
 1pc 66x43cm Non-woven cloth 1ply
 1pc Umbilical cord clamp
 2pcs 8cm Cotton sticker
 1pair PE glove
 1pair Latex examination glove
 5pcs 8x10cm Gauze swab
 1pc Face mask
 1pc Cap



G022 PE Shoe Cover

Material: PE
Color: Blue and Green
Fabric: 1.6g, 1.8g, 2.0g, 2.1g, 2.5g
 Made of PE with comfortable elastic opening.
 Available in different thickness and color.
Size: 36x15cm, 41x15cm, 40x16cm
Thickness: 0.0216mm, 0.018mm, 0.02mm, 0.021mm, 0.025mm
Packing: 100pcs/poly bag, 5000pcs/ctn
Carton Size: 48x32x36cm



G023 CPE Shoe Cover

Material: CPE
Color: Blue and Green
Fabric: 1.8g, 2.0g, 2.1g, 2.5g, 3.0g
 With comfortable elastic opening
 Available in different thickness and color
Size: 36x15cm, 41x15cm, 40x16cm
Thickness: 0.02mm, 0.025mm, 0.28mm, 0.035mm
Packing: 100pcs/Poly bag, 5000pcs/ctn
Carton Size: 48x32x36cm



G017 Lab Coat

Material: Non Woven PP
Fabric: 40g, 45g
 71x62cm Short sleeve top
 4 Button front closure
 108x56cm Pant with elastic band
Color: Blue
 Keeping the environment free of infection and
 Disease is always a battle.
 For a health professional.
 Suited to hospital and daily home use.
Packing: 1pcs/Poly bag, 50boxes/ctn
Carton Size: 56x35x50cm



G019 Overall Gown

Material: PP/PP+PE Film/SMS
 With head hood and shoe cover and zipper
Color: White, blue, Green, Dark
Weight: 35/45/50/60GSM,
Size: S, M, L, XL, XXL
 S: 164x126cm M: 172x130cm L: 178x134cm
 XL: 182x138cm XXL: 184x142cm
 Suited to lab, industrial, agriculture, hospital and daily home use.
 Unique one piece design that cover whole body convenience.
 Available in standard or PE coated.
 Coverall, zipper with storm flap attached hood, boots elastic wrists,
 waists and ankles are available bound seam construction.
Packing: 1pc/Poly bag, 50bags/ctn
Carton Size: 48x34x41cm



G025 Non-Woven Shoe Cover

Materials: PP Non-woven, SMS
Color: White and blue
Type A: PP regular shoe covers
 Comfortable elastic hem, helps you keep your cool.
 Easily worn under a sterile gown.
Type B: PP non-skid shoe covers
 Non-skid shoe covers are designed for those situations where wet or
 slippery conditions exist.
 A proprietary plastic resin is applied to the bottom of the cover to
 improve the wearers traction on the floor.
Weight: 30/40/45g
Size: 36x15cm, 38x18cm, 40x16cm
Packing: 100pcs/bag, 1000pcs/ctn
Carton Size: 52x28x350cm



G021 PE Sleeve

Material: PE
Size: 20x40cm, 42x22cm, 46x22cm
Thickness: 0.02mm/0.018mm
 Elasticized at both ends to prevent slippage fluid
 resistant.
 Effective Protection For Wrist and Forearm Clear,
 White, Blue, Red, etc.
Packing: 100pcs/poly bag, 2000pcs/ctn
Carton Size: 34x23x32cm



G029 PE Apron

Material: Polyethylene
With adjustable tie.
Available in different size and color.
Thickness: 0.018 mm, 0.020mm, 0.025mm
Size: 110x70cm, 135x71cm, 160x90cm, etc.
Color: White and blue
Packing: 100pcs/Bag, 1000pcs/ctn
Packing: 10pcs/Poly bag, 1000pcs/ctn
Carton Size: 36x25x18cm



G026 Underpad

Material: Non-Woven fabric+Fluff pulp+SAP+Fluff pulp+PE Film
Color: Blue, white, green
Size: 60x60cm, 60x90cm
SAP(Super Absorbent Polymer) Weight: 0g, 3g, 5g at customer's option
Packing: 10pcs/Bag, 120pcs/ctn
Carton Size: 50x41.5x31.5cm



G028 Non-Woven Fabrics Bedsheet

Material: Non-woven PP
Soft touch and harmless to your skin.
100% biodegradable and environmentally friendly.
SPEC: 60x180cm, 90x190cm, 100x190cm, 100x200cm, 180x200cm, 50cmx100m, other size and specification are also acceptable
Packing: 25pcs/Poly bag, 4bags/ctn
Carton Size: 32x25x24cm



G1028 Bed Sheet

Material: Film+paper
Size: 60x180cm, 90x190cm, 100x190cm, 100x200cm, etc.
Color: Blue, green, etc.
Packing: 10pcs/Poly bag, 100pcs/ctn



G027 Adult Diaper

Material:
Inner layer: Polyolefin Non-Woven
Absorbent core, dry and comfortable.
Refastening Frontal tape is convenient to put on and removed
PE film makes the surface dry and comfort
Super soft non-woven surface
Leak Guard prevent leakage efficiently
Size:
65x80cm with leg cuffs waistline 71-121cm
65x80cm without leg cuffs waistline 71-111cm
80x96cm with leg cuffs waistline 96-147cm
80x96cm without leg cuffs waistline: 96-137cm
Packing: 10pcs/Bag, 80pcs/ctn
Carton Size: 68x25.5x48cm