

H MEDICAL STERILIZATION PACKAGING



H
58

MEDICAL STERILIZATION PACKAGING

www.viormedical.org



H001 Self Sealing Flat Pouch

Size (mm)	pcs/box	pcs/ctn	Carton Size (cm)	G.W./N.W. (kgs)
57x100	200	14400	27.6x37x30	13.5/12.5
90x230	200	4000	27.3x38.2x30	14.2/13.2
135x255	200	3000	29.8x40.4x30	17/16
190x330	200	2000	37.3x39.6x30	19.7x18.7
250x370	200	1000	26.3x41.3x30	15/14



H002 Heat-Sealing Flat Pouch

Size (mm)	pcs/box	pcs/ctn	Carton Size (cm)	G.W./N.W. (kgs)
57x130	200	14400	27.6x37x30	15.8/14.8
90x230	200	4800	38.2x23x30	16/15
135x255	200	3600	26.8x42.4x30	18.1/17.1
190x330	200	2400	34.3x29.6x30	21.5x20.5
250x370	200	1200	38.3x26.3x30	16.4/15.4



H003 Heat-Sealing Gusseted Pouch

Size (mm)	pcs/box	pcs/ctn	Carton Size (cm)	G.W./N.W. (kgs)
(75+25)x300	100	3000	31.3x40x30	12.12/11.12
(100+50)x300	100	2400	31.3x42.2x30	12.8/11.8
(150+50)x400	100	1200	41.3x31.6x30	12.8/11.8
(200+55)x400	100	600	21.3x41.3x30	9.25/8.25
(250+60)x450	100	600	26.3x46.3x30	12.15/11.15



H005 Heat-Sealing Flat Reel

Size	roll/box	rolls/ctn	Carton Size (cm)	G.W./N.W. (kgs)
55mmx200m	1	6	22.5x42x23	10.58/9.58
75mmx200m	1	5	22.5x46x23	11.75/10.75
100mmx200m	1	4	22.5x47x23	10/9
150mmx200m	1	2	22.5x34x23	9.8/8.8
400mmx200m	non-box	1	22x41.5x22	12.4/11.4



H009 Autoclave Indicator Tape

Commodity: Autoclave steam tapes.

Adhesive: Acrylic pressure sensitive adhesive

Backing Material: Good quality creped paper

Specification: Size length and width are as customer request

Color Range: Light yellow to black after sterilized

Size: 12.5mmx50m, 19mmx50m, 25mmx50m, etc.



H011 Sterilization Indicator Tape For EO

Material: texture paper + EO indicators

Size: 19mmx50m, 25mmx50m

Color Range: orange to green after sterilized

Package: 1roll/PE bag



H010 Bowie & Dick Test Pack

Material: paper card

Color Change: blue to black

Size: 115mmx125mm

Package: 100pcs/ctn

Carton Size: 42x27x25cm