

D

HYPODERMIC PRODUCTS



D001 Disposable Infusion Set

Material: Medical grade PVC
 Components: Small Chamber+Plastic needle+Flow regulator+Rubber gum bubble+With latex tube+With filter
 Color: Transparent
 Sterilization: EO
 Packing: 1pc/PE Bag or blister bag, 500pcs/ctn
 Carton Size: 58x38x38cm



D1001 Disposable Infusion Set

Material: Medical grade PVC
 Components: Small chamber+Stainless needle+Flow regulator+Rubber gum bubble+With latex tube+With filter
 Color: Transparent
 Sterilization: EO
 Packing: 1pc/PE Bag or blister bag, 500pcs/ctn
 Carton Size: 58x38x38cm



D002 Disposable Infusion Set

Material: Medical grade PVC
 Components: Small chamber+Plastic needle+Flow regulator+With Y site+With filter+With luer lock adapter
 Color: Transparent
 Sterilization: EO
 Packing: 1pc/PE Bag or blister bag, 500pcs/ctn
 Carton Size: 58x38x38cm



D003 Disposable Infusion Set (Dark Color)

Material: Medical grade PVC
 Components: Small chamber+Plastic needle+Flow regulator+Rubber gum bubble+With latex tube+With filter
 Color: Brown
 Sterilization: EO
 Packing: 1pc/PE Bag or blister bag, 500pcs/ctn
 Carton Size: 58x38x38cm




D008 Disposable Infusion Set With Regulator

Material: Medical grade PVC
 Components: Small Chamber+Plastic needle+Flow regulator+With Y site+IV regulator+With Latex Tube+With Filter
 Color: Transparent
 Sterilization: EO
 Packing: 1pc/PE Bag or blister bag, 200pcs/ctn
 Carton Size: 41.5x38x45.5cm



D005 Burette Infusion Set

Material: Medical grade PVC
 Components: Air filter+Soft tube+Drip chamber+Flow regulator+Fiuld filter+Clamp+Stop valve+Y injection site+Protective cap
 Size: 100ml, 150ml, 20drops/ml, 60drops/ml
 Color: Transparent
 Sterilization: EO
 Packing: 1pc/PE Bag or blister bag, 100pcs/ctn
 Carton Size: 100ml 56.5x34.5x38cm, 150ml 62x38x38cm



D006 Disposable Blood Transfusion Set

Material: Medical grade PVC
Components: Vented spike+Drip chamber+Membrane filter+Flow regulator+Latex tube+Luer lock connector+Needle
Options available:
 With or without air vented spike,
 With or without needle,
 With or without "Y" injection port,
 Luer lock or luer slip connector.
Color: Transparent
Packing: 1pc/PE Bag or blister bag, 500pcs/ctn
Carton Size: 58x38x38cm



D1006 Disposable Blood Transfusion Set With Regulator

Material: Medical grade PVC
Components: Vented spike+Drip chamber+Membrane filter+IV regulator+Latex tube+Luer lock connector
Options available:
 With or without air vented spike; With or without needle,
 With or without "Y" injection port,
 Luer lock or luer slip connector.
Color: Transparent
Packing: 1pc/PE Bag or blister bag, 500pcs/ctn
Carton Size: 58x38x38cm



D011 Disposable Syringe

Material: Medical grade PP
2 Parts: Barrel+Plunger+Needle (or without needle)
3 Parts: Barrel+Piston+Plunger+Needle (or without needle)
Type: Luer slip, luer lock
Size: 1ml, 2ml, 3ml, 5ml, 10ml, 20ml, 30ml, 50ml, 60ml
Packing: 1pc/PE Bag or blister bag
Carton Size: PE Bag packing as below:

1ml	3600pcs/ctn	58x45x44.5cm
2ml/3ml	2700pcs/ctn	58x45x44.5cm
5ml	2400pcs/ctn	58x45x44.5cm
10ml	1600pcs/ctn	58x43x44.5cm
20ml	900pcs/ctn	58x43x44.5cm
30ml	640pcs/ctn	58x43x44.5cm
50ml/60ml	300pcs/ctn	51x43x46cm



D027 Oral Syringe

Material: Medical grade PP
Type: Luer slip, luer lock
Size: 1ml, 2ml, 3ml, 5ml, 10ml, 20ml, 35ml, 50ml, 60ml, 100ml
Package: 1pc/blister bag
Carton Size:

1ml	3600pcs/ctn	59.5x42.5x43cm
2ml/3ml	3200pcs/ctn	56.5x45x53cm
5ml	2400pcs/ctn	58x47x55cm
10ml	1500pcs/ctn	60x41x53.5cm
20ml	900pcs/ctn	56x43.5x53.5cm
50ml/60ml	200pcs/ctn	47x38.5x37.8cm
100ml	60pcs/ctn	45x29x27.5cm



D024 Disposable Infusion Pump

Material: PVC
 Infusion data
Nominal volume(ml): 60,100, 150, 200, 275
Nominal flow rate(ml/hr): 2, 3, 4, 5, 6, 7, 8, 10, 11, 12
Bolus means (ml/every time): 0.5
Lock-out time(min): 15
Package: 1pc/blister bag, 60pcs/ctn
Carton Size: 56x45x48cm



D010 Tubercle Bacillus Syringe

Material: Medical grade PP
Components: Barrel+Piston+Plunger+Needle (or without needle)
Type: Luer slip, luer lock
Size: 1ml
Packing: 1pc/PE Bag or blister bag, 3600pcs/ctn
Carton Size: 16.5x48.5x38cm



D1009 Irrigation Syringe

Material: Medical grade PP
Components: Barrel+Piston+Plunger
Type: Cap or without cap
Size: 50ml, 60ml, 100ml, 150ml, 200ml
Sterilization: EO
Packing: 1pc/PE Bag or blister bag
Carton Size: 50/60ml 300pcs/ctn 51x43x46cm



D1012 Irrigation Syringe With Ball

Material: Medical grade PVC/PP material
Color: Blue and transparent
Size: 60cc (ml)
Sterilization: EO
Package: 1pc/Blister bag, 50pcs/ctn
Carton Size: 48x38x23cm



D009 Auto Destruct Syringe

Material: Medical grade PP
Components: Barrel+Piston+Plunger+Clip+Fixed Needle
Size: 0.2ml, 0.3ml, 0.5ml, 1ml
Packing: 1pc/blister bag, 100pcs/box, 1800pcs/carton
Carton Size: 48.5x39.5x28.5cm



D034 Auto Disable Syringe

Material: Medical grade PP
 Back-end structure
Components: Barrel+Piston+Plunger+Needle (or without needle)
Size: 2ml, 3ml, 5ml, 10ml
Packing: 1pc/PE Bag or blister bag
Carton Size:

2ml/3ml	2400pcs/ctn	58x44x46.5cm
5ml	1800pcs/ctn	58x44x46.5cm
10ml	1200pcs/ctn	56x44x46cm



D012 Insulin Syringe

Material: Medical grade PP
Components: Barrel+Piston+Plunger+Needle+Cap
Type: 40IU, 100IU
Size: 0.3ml, 0.5ml, 1ml
Needle Size: 29G, 30G, 31G
Packing: 1pc/PE Bag or blister bag
Carton Size:

0.3ml/0.5ml	200pcs/Box, 3000pcs/ctn	48.5x40.5x43.5cm
1ml	100pcs/Box, 3200pcs/ctn	61.5x48.5x38cm



D020 Insulin Pen Needle

Material: PP, Stainless steel
Needle size: 29G, 30G, 31G, 32G, 33G
Needle length: 4mm, 6mm, 8mm
Sterilization: EO
Package: 1pc/hard blister bag, 100pcs/box, 5000pcs/ctn
Carton Size: 37.5x27x27.5cm



D013 Disposable Needle

Material: Stainless steel
Components: PP base+Needle
Size: 18G, 19G, 20G, 21G, 22G, 23G, 24G, 25G, 26G, 27G, 28G, 29G, 30G
Color: Pink, Cream, Yellow, Green, Black, Blue, Purple, Orange, Brown, Gray, Blue-green, Red, Light-yellow
Packing: 1pc/Blister bag, 100pcs/Box, 100box/ctn
Carton Size: 45x46x40cm




D033 Disposable Safety Needle

Material: Stainless steel
Size: 18G, 19G, 20G, 21G, 22G, 23G, 24G, 25G, 26G, 27G, 28G, 29G, 30G
Color: Pink, Cream, Yellow, Green, Black, Blue, Purple, Orange, Brown, Gray, Blue-green, Red, Light-yellow
Packing: 1pc/Blister bag, 100pcs/Box, 2000pcs/ctn
Carton Size: 50.5x43.5x40.5cm



D025 Epidural Spinal Needle Set

Components: Epidural needle(tuohy needle), 16G/18Gx80mm
Spinal needle: 25Gx110mm, pencil point
Epidural catheter: multiport with mark 0.8/1.0mm850mm
Epidural filter: 0.22m, hydrophobic membrane, Luer lock
LOR syringe, 7ml, 10ml available
Catheter adaptor: TBA, secure the catheter completely.
Packing: 1set/Hard Blister bag, 100sets/ctn
Carton Size: 56x34x29cm



D021 Central Venous Catheter Kit(CVC Kit)

Material: TPU, silicone, stainless steel
Model Size: 14G, 16G, 18G, 20G
Type: Single Lumen, double Lumen, Triple Lumen
Model Size: 4F, 5F, 7F, 8F
Packing: 1set/Hard Blister bag, 10 Sets/Box, 120sets/ctn
Carton Size: 48x42x76cm



D2009 Comestic Syringe

Material: Medical grade PP
Components: Barrel+Piston+Plunger
Color: Transparent
Size: 1ml, 3ml
Sterilization: EO
Packing: 1pc/PE Bag or blister bag
Carton Size:
 1ml 3600pcs/ctn 58x45x44.5cm
 3ml 2700pcs/ctn 58x45x44.5cm



D028 Fine Micro Cannula

Material: PP, Stainless Steel
 Blunt tip for a non-traumatic introduction and reducing risk of hematoma
 Color coded according to the diameter of the cannula
 Lateral opening on the tip ensures a precise injection
 With graduation or without graduation
Size: 18G, 21G, 23G, 25G, 27G, 30G
Needle length: 38mm, 50mm, 70mm
Sterilization: EO
Package: 1pc/Blister bag, 50pcs/Box, 2500pcs/ctn
Carton Size: 46x42x49.5cm



D007 Stopcock Three Way

Material: PC+PE+ABS
Components: Body+Plug+Screw cap
Options available: With tube or without tube
Products Thread: 6:100 screw thread standard
 This in line with international GB standards
Packing: 1pc/Hard blister, 100pcs/box, 1000pcs/ctn
Carton Size: 62x32x37cm



D022 Extension Tube

Material: PVC, PP
Size: 150cmx3.0mm, 30cmx3.0mm, 60cmx3.0mm
Package: 1pc/blister bag, 50pcs/Box, 500pcs/ctn
Carton Size: 77x32x34cm



D015 Spinal Needle

Material: Stainless steel
Components: PP hub+Needle
Type: Pencil point tip, Quincke tip
Size: 18G, 19G, 20G, 21G, 22G, 23G, 24G, 25G, 26G, 27G
Color: Pink, Cream, Yellow, Green, Black, Blue, Purple, Orange, Brown, Gray
Length: 90mm-220mm
Packing: 1pc/Blister bag, 1000pcs/ctn
Carton Size: 42x27x52cm



D016 Epidural Needle

Material: Stainless steel
Components: PP hub+Needle
Size: 15-18G
Length: 80mm, 90 mm
Packing: 1pc/Blister bag, 50pcs/Box, 1000pcs/ctn
Carton Size: 54x42x36cm



D023 Combi Stopper

Material: PVC, PP, ABS
Sterilization: EO
Color: Red, blue
Package: 1pc/blister bag, 100pcs/box, 5000pcs/ctn
Carton Size: 42X26X45cm



D029 Needle Free Connector

Material: PC
Sterilization: EO
Package: 1pc/hard blister bag, 50pcs/box, 2500pcs/ctn
Carton Size: 42x26x45cm



D1004 I.V. Flow Regulator

Material: PVC, PP, ABS
Type: Single graduation, Double graduation
 With tube or without tube
Sterilization: EO
Package: 20pcs/box, 400pcs/ctn
Carton Size: 40x40x30cm



D030 Heparin Cap

Material: PC+isoprene rubber, ABS+isoprene rubber
Sterilization: EO
Package: 1pc/hard blister bag, 100pcs/box, 5000pcs/ctn
Carton Size: 42x26x45cm



D017 Scalp Vein Set

Material: Medical grade PVC+Stainless steel
Wing Type: Butterfly wing, Single wing
Size: 18G, 19G, 20G, 21G, 22G, 23G, 24G, 25G, 26G, 27G
Color: Pink, Cream, Yellow, Green, Black, Blue, Purple, Orange, Brown, Gray
Length: 15-30cm
Packing: 1pc/PE Bag or blister bag, 100pcs/box, 2000pcs/ctn
Carton Size: 55x37x26cm



D1010 I.V. Cannula

Material: Stainless steel
Type: With butterfly wings, Injection port with wings, Pen-like, Safety Pen-like
Size: 14G, 16G, 18G, 20G, 22G, 24G, 26G
Packing: 1pc/Blister bag,
Carton Size:
Injection port: 1000pcs/ctn 72x31x31.5cm
Pen-like: 1000pcs/ctn 47.5x38x25cm



D019 A.V. Needle

Material: Stainless steel
Size: 15G, 16G, 17G
Length: 30cm
Packing: 1pc/PE Bag or blister bag, 50pcs/box, 500pcs/ctn
Carton Size: 63x36x45.5cm



D1011 Safety I.V. Cannula

Size: 18G, 20G, 22G, 24G, 26G
Color: Green, Pink, Blue, Yellow, Purple
Packing: 1pc/Hard blister, 200pcs/box, 1200pcs/ctn
Carton Size: 58x42x47cm



D032 Acupuncture Needle

Material: 304 Medical Stainless Steel
Handle: Stainless Steel Wire Handle
Needle diameter: 0.12-0.35mm
Length: 13mm-100mm
Tube: one needle with one tube
Package: Sleeve Packaging, 5pcs/hard blister bag, 100pcs/box, 10000pcs/ctn
Carton Size: 60x31x25cm



D3009 Blood Bag

Type: Extruding bag, Rolled bag
Material: Medical grade PVC
Color: White transparent
Packing: 1pc/PE Bag

Specification Bag	Capacity	Anticoagulant	QTY	Carton Size
Single	250ml-500ml	CPDA/CPD	100	51x32x26cm
Double	250ml-500ml	CPDA/CPD	60	51x32x26cm
Triple	250ml-500ml	CPDA/CPD+SAG-M	32	51x32x24cm
Quadruple	250ml-500ml	CPDA/CPD+SAG-M	30	36x36x20cm
Quintuple	250ml-500ml	CPDA/CPD+SAG-M	18	51x32x20cm
Transfer bag	150ml, 300ml	Empty	100	51x32x26cm