

B MEDICAL TUBES

B
11

B013 Closed Suction Catheter

Material: PVC, PP, ABS

Application: Used with Endotracheal tube or Tracheotomy tube to clear bacterial buildup and reduce VAP pathogenesis.

Type: 24h, 72h

Size: Fr6, Fr8, Fr10, Fr12, Fr14, Fr16

Color coded control valve for easily identify sizes.

Sterilization: EO

Packing: 1pc/ Blister bag, 200pcs/ctn

Carton Size: 74x39x36cm



B001 Suction Catheter

Material: Non-toxic, non-irritant soft PVC

Type of connector:

T-type connector, Y-type connector,

Cap cone connector, Plain type connector

With X-ray or without is available

Standard Tube Length: 50cm

Size: Fr5-Fr24

Packing: 1pc/PE Bag, 600pcs/ctn

Carton Size: 57x37x27cm

Packing: 1pc/Blister bag, 600pcs/ctn

Carton Size: 57x42x31cm



B002 Extractor Mucus

Material: Soft, frosted and kink resistant PVC tubing

Spare screw top lid is provided to seal the

container for safe transportation of specimen

to the laboratory or aseptic disposal of container

Suction tube length available: 40cm, 50cm

Capacity: 25ml, 40ml

Size: Fr6-Fr14

Packing: 1pc/PE Bag, 300pcs/ctn

Carton Size: 44x42x45cm



B005 Stomach Tube

Material: Non-toxic, non-irritant soft

PVC or made from 100% silicone

Printed with graduation.

Color coded for size identification

With X-ray or without is available

Size: Fr6-Fr24

Tube Length: 120cm

Packing: 1pc/PE Bag

Carton Size: 57x37x27cm

Packing: 1pc/Blister bag

Carton size: 50x42x31cm

Size	QTY	Size	QTY
Fr6 - Fr8	600	Fr12 - Fr14	360
Fr10	480	Fr16 - Fr18	300



B014 Latex Suction Catheter

Material: Made of Natural Latex material, silicone coated


Type of connector: T-type connector, Y-type connector,

Cap cone connector, Plain type connector

Size: 6Fr-22Fr

Package: 1pc/blister bag, 50 pcs/box, 20boxes/ctn

Carton Size: 73x52x36cm



B003 Feeding Tube

Material: Non-toxic, non-irritant soft PVC

Color coded for size identification

With X-ray or without is available

Size: Fr4-Fr8

Tube Length: 40cm

Size: Fr10-Fr22

Tube Length: 120cm

Packing: 1pc/PE Bag

Carton Size: 57x37x27cm

Packing: 1pc/Blister bag

Carton size: 50x42x31cm

Size	QTY	Size	QTY
Fr4 - Fr8	1200	Fr12 - Fr14	360
Fr10	480	Fr16 - Fr18	300



B016 Silicone Stomach Tube

Material: Silicone

Made of 100% medical grade silicone

X-ray opaque line throughout the length

Open distal end with 4 side holes

Size: Fr10-Fr20

Package: 1pc/Blister Bag, 20pcs/Box, 160pcs/ctn

Carton Size: 60x41x35cm

MEDICAL TUBES

VIOR
MEDICAL

B
12

MEDICAL TUBES



B0070 Nelaton Catheter(PVC)

Material: Non-toxic, non-irritant soft PVC
Color coded for size identification
With X-ray or without is available
Size: Fr6-Fr24 Tube Length: 40cm
Packing: PE Bag, 600pcs/ctn Carton Size: 57x37x27cm
Packing: Blister bag, 600pcs/ctn Carton Size: 50x42x31cm



B006 Rectal Tube

Material: Non-toxic, non-irritant soft PVC
Color coded for size identification
With X-ray or without is available
Size: Fr20-Fr30 Tube Length: 40cm
Packing: PE Bag, 500pcs/ctn
Carton Size: 57x37x27cm



B008 Penrose Drainage Tube

Material: Natural latex
Size: 1/4", 3/8", 1/2", 5/8", 3/4", 7/8", 1"
Length: 12", 18"
Packing: 1pc/Blister bag, 1000pcs/ctn
Carton Size: 47x25x54cm



B009 T-Drainage Tube

Material: Natural latex, silicone coated
Used for the drainage of the operation of Cholecyst.
The connector part is larger in diameter than other parts.
Size: Fr8-Fr28
Packing: 1pc/PE Bag, 200pcs/ctn
Carton Size: 49.5x37x34cm



B0072 Lubricated Nelaton Catheter

Material: Non-toxic, non-irritant soft PVC or PU
Color coded for size identification. With X-ray or without is available
Size: Fr6-Fr18

Type	French Size	Length	Tip type	PCS/Box
Male	8Fr, 10Fr, 12Fr, 14Fr, 16Fr, 18Fr	40cm	Straight	30pcs
Coude Male	10Fr, 12Fr, 14Fr, 16Fr	40cm	Tapered or Olive	30pcs
Pediatric	6Fr, 8Fr, 10Fr, 12Fr	25cm	Straight	30pcs
Female	8Fr, 10Fr, 12Fr, 14Fr, 16Fr	19cm	Straight	30pcs



B010 Umbilical Catheter

Material: Medical grade PVC.
Designed for intermittent continuance access to the umbilical artery or vein of newly born or premature infant.
Size: FG-5(Gray), FG-6(Green) and FG-8(Blue)
Length: 125cm
Tube is marked at 5, 6, 7, 8, 9, 10, 15cm
Packing: 1pc/PE Bag, 125pcs/box, 2500pcs/ctn
Carton Size: 44x42x46cm



B011 Condom Catheter

Material: Natural latex
Size: 20mm, 25mm, 30mm, 35mm, 40mm
Packing: 1pc/Blister bag, 100pcs/box, 1000pcs/ctn
Carton Size: 56x54x28cm



B012 Condom Catheter

Material: 100% silicone medical grade
Size: 18mm, 24mm, 28mm, 31mm, 35mm, 40mm
With/Without adhesive bag
Packing: 1pc/PE Bag, 30pcs/box, 300pcs/ctn
Carton Size: 78x39.5x28cm