




E

SURGERY & GYNECOLOGICAL EXAMINATION PRODUCTS



E

30

SURGERY & GYNECOLOGICAL
EXAMINATION PRODUCTS

www.viormedical.org



E001 Surgical Gloves

Material: Natural Latex
Type: Powdered or Powder-free
Sterilization: Gamma
Packing: 50pairs/box, 500pairs/ctn
Carton Size: 69x25x40cm

Size	6	6.5	7	7.5	8	8.5	9
Length(min, mm)	260	260	270	270	270	280	280
Width(min, mm)	77±5	83±5	89±5	95±5	105±6	108±6	114±5
(70°Cx168hours) tensile Strength(mpa)	Before Aging			After Aging			
	≥12			9			



E003 Nitrile Examination Gloves

Material: Nitrile
Size: XL, L, M, S
Type: Powdered or Powder-free
Packing: 100pcs/box, 1000pcs/ctn
Carton Size: 30.5x26x24cm

Size (±10mm)	XL	L	M	S
Length(min, mm)	230			
Width(min, mm)	110	100	95	80
Thickness Single (mm)	≥0.08			
Physical Properties(mm)	26pcs Sampled			
(100°Cx22hours)tensile Strength(mpa)	Before Aging		After Aging	
	14.0		13.0	
Ultimate Elongation(%)	520		410	



E004 Latex Gynaecological Gloves

Material: Natural Latex
Type: Powdered or Powder-free
Size: L, M, S
Sterilization: Gamma
Packing: 25pairs/box, 100pairs/ctn
Carton Size: 38x36x38cm

Dimensions		Standards	
		Yingmed Gloves	EN455
Length(mm)		450	400
Width(mm)	Small	83±6	83±6
	Medium	95±6	95±6
	Large	108±6	108±6
Thickness Single wall(mm)	Fingers	>0.2	>0.1
	Palm	>0.16	>0.1
	Cuff	>0.12	>0.1



E005 Vinyl Gloves

Material: Vinyl
Grade: AQL2.5
Length: 240±10mm
Thickness: 0.08±0.01mm
Weight: 4.5g±0.2
Content of Powderful: 0.1-0.2g
Size: XL, L, M, S
Width Of The Centre Of The Paim:
XL 11.4±0.3cm, M 10.8cm±0.3, S 10.6±0.3cm
Packing: 100pcs/box, 1000pcs/ctn 31.5x23x22cm

YM-E006 PE Gloves

Material: Recycled or virgin polyethylene
Thickness: 0.018mm, 0.02mm, 0.022mm etc.
Size: L, M, S
Type: HDPE or LDPE, plain or embossed
Color: White, blue, pink, etc.
Packing: 100pcs/Box or bag, 10000pcs/ctn
Carton Size: 40x32x26cm
Different thickness, sizes, color and packings is available.

YM-E007 Two Fingers Cots

Weight: 1.2g/pc
Color: white
Packing: 100pcs/bag, 10000pcs/ctn
Carton Size: 39x30x24cm



E002 Latex Examination Gloves

Material: Natural Latex
Type: Powdered or Powder-free
Size: L, M, S
Packing: 100pcs/box,1000pcs/ctn
Carton Size: 37x25x24cm

Size	L	M	S	±10mm
Length(min, mm)	230			
Width(min, mm)	110	95	80	/
Thickness Single (mm)	≥0.08			
Physical Properties(mm)	26pcs Sampled			
(100°Cx22hours)tensile Strength(mpa)	Before Aging		After Aging	
	21.0		16.0	
Ultimate Elongation(%)	700		500	





E008 Disposable Surgical Blade

Material: Carbon steel or stainless steel
Surgical Blade are critically made in accordance with ISO9001/ISO7740 standard. They supply very high quality to high level doctors world wide.
Size: 10#-36#
Packing: 1pc/Aluminium foil bag, 100pcs/box, 5000pcs/ctn
Carton Size: 36x32x10.5cm



E009 Disposable Skin Graft Blade

Material: Carbon steel or stainless steel
Size: 160mm
Packing: 1pc/blister bag, 2500pcs/ctn
Carton Size: 43.5x42.5x23.5cm



E013 Stainless Steel Lancet

Material: Stainless steel
Size:
L(40x5.5x0.15mm), S(32x5x0.15mm)
Packing: 1pc/paper bag, 200pcs/box, 100boxes/ctn
Carton Size: 44x33.5x18cm



E015 Twist Blood Lancet

Material: Stainless steel+Plastic handle
Size: 21G, 23G, 26G, 28G, 30G
Sterilization: Gamma
Packing: 100pcs/box, 200boxes/ctn
Carton Size: 46x35x20cm



E012 Disposable Surgical Scalpel

We supply many different handles for different operation requirements. Protective tab does best safety job to operators.
Disposable Scalpels(Inserted model)
Packing: 1pc/Aluminium foil bag, 10pcs/Box, 500pcs/ctn
Carton Size: 40x32x23cm



E045 Skin Stapler

Material: Metal, Plastic, Steel
Anastomotic nail number: 15W, 25W, 35W
Anastomotic nail type: R-type, W-type
Package: 1pc/blister bag, 5pcs/box, 100boxes/ctn
Carton Size: 60.5x48x21cm



E040 Safety Blood Lancet

Material: ABS+Stainless steel
Pressure type: Gamma sterile
Size: 18G, 21G, 23G, 26G, 30G
Packing: 100pcs/Box, 20box/ctn
Carton Size: 33x28.4x15.5cm



E019 Lancing Device

Lancing device with ejector
Depth: 5
Lancing Device assist diabete patients to monitor their glucose levels, easily and safely with an adjustable dial, lancing depth can be set to an appropriate level for the individual. More safe professional system. Any standard lancet can be used with these devices.
Packing: 1pc/box, 100pcs/ctn
Carton Size: 27.5x22x22cm



E044 Surgical Skin Marker

Material: Non Toxic PP
Size: 0.5mm, 1mm, double tip or single tip
Ink: Dry-Erase & Wet-Erase, white/purple/blue
Specification: 14cm pen with 15cm soft plastic ruler
Sterilization: Gamma
Package: 1pc/blister bag, 100pcs/box, 800pcs/Carton
Carton Size: 49x35x35cm



E046 Stapler Remover

Material: Metal, Plastic
Type: for prepuce, other derma
Package: 1pc/blister bag



E042 Conductive Gel

Material: Aqueous matrix polymer gel with electrolyte
PH value: 5.5-9.0
Viscosity: 15000~80000mPa.S
Electrical conductivity(25°C): 10.0~85.0ms/cm
Capacity: 100ml, 250ml, 500ml
Package: 100ml/bottl, 100bottles/40x25x23cm



E043 Ultrasound Transmission Gel

Size: 250ML, 1L, 2.5L
Color: Clear or light blue
Package: 250ml/bottl, 40bottles/41x26x20cm
5L/bottl, 2bottles/38X20X24cm



E018 Disposable Medical Razor

Carbon steel or stainless steel, single or double blade.
All razors come with the highest quality Stainless Steel
Packing: 100pcs/box, 1000pcs/ctn
Carton Size: 63x32x21cm



E047 Lubricating Jelly

Material: Gel water base, for medical use
Size: 3g, 5g, 10g, 30g, 42g, 45g, 82g, 113g, ect
Package: 1pc/Box, 72pcs/ctn
Carton Size: 38.6x34.5x20cm(113g)

Surgical Suture

Features:

- 1. Thread Length: 45cm, 60cm, 75cm, 100cm, 125cm, 150cm ect.
- 2. Thread diameter: 8/0, 7/0, 6/0, 5/0, 4/0, 3/0, 2/0, 1/0, 0, 1, 2, 3, ect.
- 3. Needle Curvature: straight, 1/2 circle, 1/2 circle (double), 1/4 circle, 1/4 circle (double), 3/8 circle, 3/8 circle(double), 5/8 circle, loop round.
- 4. Needle Length: 6mm, 8mm, 12mm, 18mm, 22mm, 30mm, 35mm, 40mm, 50mm etc.
- 5. Needle End Type: rolled end, drilled end, regular-eye end.
- 6. Needle Body: taper point, blunt taper point, cutting edge, reverse cutting edge, tapercut, micro-point spatula curved.
- 7. Cross Section: round bodied, round bodied(heavy), curved cutting, curved cutting(heavy).
- 8. Reverse Cutting: reverse cutting(heavy), tapercut, micro-point spatula curved etc.
- 9. Sterilization Method: gamma radiation or ethylene oxide(can be used directly without being re-sterilized before use).

Packing:

- 1 set/sealed polyester and aluminum foil sachet,
- 12 foil sachets/Box or plastic container, 50boxes/ctn 28x28x34cm



Absorbable Surgical Suture	Natural	E020	Chromic Catgut
		E021	Plain Catgut
	Artificial chemical synthesis	E022	Polyglycolic Acid
		E032	Rapid Polyglycolic Acid
		E028	Polyglactin 910 (PGLA 910)
		E038	Rapid Polyglactin 910 (PGLA 910-R)
		E029	Polydioxanone (PDO/ PDS)
Non-absorbable Surgical Suture	Natural	E030	Poleglecaprone(PGCL)
	Artificial chemical synthesis	E023	Silk Braided
		E025	Nylon Monofilament
		E026	Polyester Braided
		E027	Polypropylene
		E031	Stainless Steel Wire



Vaginal Speculum

Description:

It can be used for the checking and treatment of gynecology disease in hospitals or clinics. Before using it, user must read the usage manual, and then take the speculum, push the "duckbill" of speculum that is closed into the vaginal slowly .
Open the "duckbill" of speculum according to demand, fix the nut of speculum, then it can be used for checking and treatment of gynecology disease. After all, Loose the nut of speculum, and take "duckbill" of speculum out of the vaginal slowly.

Specification:

Material: Polystyrene, transparent substance, non-toxic, medical grade.





Size: Large, Medium, Small

Sterilization: EO

Packing: Steriled in individual polybag or Film+Polybag or Blister pack, 100pcs/ctn
Carton Size: 520x380x360mm G.W: 7kgs N.W: 6kgs

E1013 Disposable Vaginal Speculum With Side Screw
E1016 Disposable Vaginal Speculum Spanish Type
E1018 Disposable Vaginal Speculum Australia Type

E1015 Disposable Vaginal Speculum French Type
E1017 Disposable Vaginal Speculum With Middle Screw
E1019 Disposable Vaginal Speculum South Asian Type



E1021 Cervical Rambrush

Usage: Sampling Examination For Male And Female
Material: Stainless Steel Thread + PVC Handle + Nylon Tip
PVC Model: Slid Type, Injection Type
Packing: 1pc/blister bag
Carton Size: Normal Type 5000pcs/ctn 41x35x14cm
Broom Type 4000pcs/ctn 44x43x37cm



E1022 Amniotomy Hook

Material: Medical grade PVC
Packing: 1pc/blister bag, 100pcs/box, 20box/ctn
Carton Size: 52x50x30cm



E1023 Cervical Scraper

Material: wooden
Packing: 100pcs/box, 50box/ctn
Carton Size: 46.5x37.5x17cm



E1024 Cervical Scraper

Material: plastic
Packing: 2500pcs/ctn
Carton Size: 57x25x50cm



E2009 Hcg Test Devise

Size: 3.0mm/3.5mm/4.0mm
Specimen: Urine, Urine/Serum
Sensitivity: 10mIU/ml, 15mIU/ml, 20mIU/ml, 25mIU/ml
Packing: 2500pcs/ctn 75x44x39cm



E2010 Hcg Midstream

Size: 3.0mm/4.0mm/5.5mm/6.0mm/7.0mm
Specimen: Urine, Urine/Serum
Sensitivity: 10mIU/ml, 15mIU/ml, 20mIU/ml, 25mIU/ml
Packing: 2000pcs/ctn 75x39x44cm



E2017 HCG Test Strip

Size: 2.5mm
Packing: 1pc/foil pouch, 50pcs/box, 5000pcs/ctn
Carton Size: 75x29x56cm



E1025 Disposable PAP Smear Kit

Content: 1pc X Wooden Cervical scraper, 1pc X Cervical brush
1pc X Slide glass, 1pc X Cotton tipped applicators
Packing: 1set/paper box/Blister bag, 400sets/ ctn 57.5x43x22cm