

C

UROLOGY PRODUCT



C
14

UROLOGY PRODUCT

C001 2-Way Latex Foley Catheter

Material: Natural latex with silicone coated

Type: Plastic valve, Rubber valve

Specifications:

Fr6-Fr10 (Pediatric)	Length: 270mm	Balloon Capacity: 3cc/5cc
Fr12-Fr24(Standard)	Length: 400mm	Balloon Capacity: 5cc/10cc/30cc
Fr26-Fr30(Standard)	Length: 400mm	Balloon Capacity: 10cc/30cc

Color: Coded for visualization of size

Packing: 1pc/Blister Bag, 10pcs/box, 500pcs/ctn


Carton Size: 52x42x47.5cm

www.viormedical.org



C002 3-Way Latex Foley Catheter

Material: Natural latex with silicone coated
 Type: Plastic valve, Rubber valve
 Specifications:
 Length: 400mm
 Fr16-Fr24(Standard) Balloon Capacity: 30cc
 Fr26-Fr28(Standard) Balloon Capacity: 10cc/30cc
 Color: Coded for visualization of size
 Packing: 1pc/Blister bag, 10pcs/box, 500pcs/ctn
 Carton Size: 52x42x47.5cm




C003 2-Way Female Foley Catheter

Material: Natural latex with silicone coated
 Type: Plastic valve, Rubber valve
 Specifications:
 Length: 240mm
 Fr12-Fr22 Balloon Capacity: 5cc/10cc/30cc
 Fr24-Fr26 Balloon Capacity: 10cc/30cc
 Color: Coded for visualization of size
 Packing: 1pc/Blister bag, 10pcs/box, 500pcs/ctn
 Carton Size: 52x42x47.5cm



C005 2-Way Tiemann Foley Catheter

Material: Natural latex with silicone coated
 Type: Plastic valve, Rubber valve
 Specifications:
 Length: 400mm
 Fr12-Fr22 Balloon Capacity: 5cc/10cc/30cc
 Fr24 Balloon Capacity: 10cc/30cc
 Color: Coded for visualization of size
 Packing: 1pc/Blister bag, 10pcs/box, 500pcs/ctn
 Carton Size: 52x42x47.5cm



C006 3-Way Double Balloons Foley Catheter

Material: Natural latex with silicone coated
 Type: Plastic valve, rubber valve
 Specifications:
 Length: 400mm
 Size: Fr16-Fr24
 Balloon Capacity: 30cc(front)/50cc(rear)
 Color: Coded for visualization of size
 Packing: 1pc/Blister bag, 10pcs/box, 500pcs/ctn
 Carton Size: 52x42x47.5cm



C008 4-Way Double Balloons Foley Catheter

Material: Natural latex with silicone coated
 Specifications:
 Size: Fr16-Fr24 Length: 400mm
 Balloon Capacity: 30cc(front)/50cc(rear)
 Color: Coded for visualization of size
 Packing: 1pc/Blister bag, 10pcs/box, 500pcs/ctn
 Carton Size: 52x42x47.5cm



C009 Urethral Catheter(Red Latex)

Material: Natural latex with silicone coated
 Specifications: Fr12-Fr40
 Length: 400mm
 Packing: 1pc/Blister bag, 10pcs/box, 500pcs/ctn
 Carton Size: 52x43x47.5cm



C010 Malecot Catheter

Material: Natural latex with silicone coated
 Specifications: Fr6-Fr24
 Length: 400mm
 Packing: 1pc/Blister bag, 10pcs/box, 500pcs/ctn
 Carton Size: 52x42x47.5cm




C060 2-Way Tiemann Silicone Foley Catheter

Material: High quality 100% silicone, medical grade
 Type: Hard valve, with X-ray or without is available
 Specifications:
 Size: Fr14-Fr26
 Color: Coded for visualization of size.
 Packing: 1pc/Blister bag, 10pcs/box, 500pcs/ctn
 Carton Size: 52x43x39.5cm



C007 2-Way All Silicone Foley Catheter

Material: high quality 100% silicone, medical grade
 Type: Hard valve, with X-ray or without is available
 Specifications:
 Color: Coded for visualization of size.
 Fr6-Fr10(Pediatric) Length: 310mm Balloon Capacity: 3cc/5cc
 Fr12-Fr18(Standard) Length: 400mm Balloon Capacity: 5cc/10cc
 Fr20-Fr24(Standard) Length: 400mm Balloon Capacity: 5cc/10cc/30cc
 Packing: 1pc/Blister bag, 10pcs/box, 500pcs/ctn
 Carton Size: 52x43x39.5cm



C061 3-Way All Silicone Foley Catheter

Material: High quality 100% silicone, medical grade
 Type: Hard valve, with X-ray or without is available
 Specifications:
 Size: Fr18-Fr24
 Color: Coded for visualization of size.
 Packing: 1pc/Blister bag, 10pcs/box, 500pcs/ctn
 Carton Size: 52x43x39.5cm



C011 Economic Urinary Drainage Bag Without Outlet

Capacity: 2000ml
Without bottom outlet, with 90cm inlet tube,
With non-return valve.
Sterilization: EO
Packing: 1pc/PE Bag, 10pcs/bag, 250pcs/ctn
Carton Size: 40x36x37cm



C012 Economic Urinary Drainage Bag With Push Pull Valve

Capacity: 2000ml
With bottom outlet, with Push Pull Valve,
With 90cm inlet tube, with non-return valve.
Sterilization: EO
Packing: 1pc/PE Bag, 250pcs/ctn
Carton Size: 40x36x37cm



C018 Luxury Urinary Drainage Bag

Capacity: 2000ml
Type: Square bottom, latex free
With T valve,
With white sample port,
With anti-reflux valve, with plastic tie
Inlet tube: O.D 10 mm, 100cm length
Sterilization: EO
Packing: 1pc/Blister pack, 100pcs/ctn
Carton Size: 54x36x32cm



C020 Luxury Urinary Drainage Bag

Capacity: 2000ml
Type: Square bottom, latex free
With Single reversal valve, with needle sample port
With anti-reflux valve, with double plastic hanger and
single rope, with blue bed sheet clip.
Inlet tube: O.D 10 mm, 100cm length
Sterilization: EO
Packing: 1pc/Blister pack, 40pcs/ctn
Carton Size: 43x32x38cm



C013 Economic Urinary Drainage Bag With T valve

Capacity: 2000ml
With bottom outlet, with T valve,
With 90cm inlet tube, with non-return valve.
Sterilization: EO
Packing: 1pc/PE Bag, 250pcs/ctn
Carton Size: 40x36x37cm



C062 Economic Urinary Drainage Bag Push Pull Screw Valve

Capacity: 2000ml
With bottom outlet, with screw valve,
With 90cm inlet tube, with non-return valve.
Sterilization: EO
Packing: 1pc/PE Bag, 250pcs/ctn
Carton Size: 40x36x37cm



C022 Urinary Drainage Bag

Capacity: 2000ml
Type: Pyriform, latex free
With T valve, with needle sample port, with anti-reflux valve,
With double plastic hanger and rope,
With white bed sheet clamp.
Inlet tube: O.D 10 mm, 110cm length
Sterilization: EO
Packing: 1pc/Blister pack, 40pcs/ctn
Carton Size: 43x30x38cm



C023 Urinary Drainage Bag

Capacity: 2000ml
Type: Pyriform, latex free
With T valve, with needle sample port, with anti-reflux valve,
with blue plastic hanger and rope,
with white bed sheet clamp.
Inlet tube: O.D 10 mm, 100cm length
Sterilization: EO
Packing: 1pc/Blister pack, 50pcs/ctn
Carton Size: 55x40x27cm



C015 Luxury Urinary Drainage Bag

Capacity: 2000ml
Type: Square bottom, latex free, with pull push valve,
with white sample port, with anti-reflux valve, with
plastic tie
Inlet tube: O.D 10 mm, 100cm length
Sterilization: EO
Packing: 1 pc/PE or Blister pack, 100pcs/ctn
Carton Size: 54x36x32cm



C017 Luxury Urinary Drainage Bag

Capacity: 2000ml
Type: Square bottom, latex free
With T valve, with white sample port, With anti-reflux valve,
without plastic tie
Inlet tube: O.D 10 mm, 100cm length
Sterilization: EO
Packing: 1pc/Blister pack, 100pcs/ctn
Carton Size: 54x36x32cm



C025 Urinary Drainage Bag

Capacity: 2000ml
Type: Pyriform, latex free
With T valve, with needle sample port,
With anti-reflux valve,
With white double plastic hanger and rope,
With white bed sheet clamp.
Inlet tube: O.D 10mm, 110cm length
Sterilization: EO
Packing: 1pc/Blister pack, 40pcs/ctn
Carton Size: 43x30x38cm



C026 Luxury Urinary Drainage Bag

Capacity: 2000ml
With T valve, Latex free,
With anti-reflux drip chamber(one part), with needle sample
port and tubing clip is optional, with anti-reflux valve.
With reinforced double plastic hanger and rope.
With blue bed sheet clamp.
Inlet tube: O.D 11mm, 110cm length
Sterilization: EO
Packing: 1pc/Blister pack, 40pcs/ctn
Carton Size: 43x30x38cm



C028 Luxury Urinary Drainage Bag

Capacity: 2000ml
 With T valve, latex free
 With anti-reflux drip chamber(three part)
 With needle sample port and tubing clip is optional
 With anti-reflux valve
 With reinforced double plastic hanger and rope
 With blue bed sheet clamp
Inlet tube: O.D 11 mm, 110cm length
Sterilization: EO
Packing: 1pc/Blister pack, 40pcs/ctn
Carton Size: 43x30x38cm



C029 Luxury Urinary Drainage Bag

Capacity: 2000ml
 With T valve, latex free
 With anti-reflux drip chamber(three part)
 With needle sample port and tubing clip is optional
 With anti-reflux valve
 With reinforced double plastic hanger and rope
 With blue bed sheet clamp
Inlet tube: O.D 11 mm, 110cm length
Sterilization: EO
Packing: 1pc/Blister pack, 40pcs/ctn
Carton Size: 43x30x38cm



C037 Disposable Pediatric Urine Collector

Size: 100ml, 200ml
Component: Pectic paper, PVC bag, sponge
Sterilization: EO
Packing: 1pc/PE Bag, 100pcs/box, 5000pcs/ctn
Carton Size: 55x54x43cm



C038 Wound Drainage System

Material: PVC and/or silicone rubber, PP, PS, SS
Application: Used for abdominal, chest, breast and other parts of the fluid, pus and blood drainage.
Component: Container, drainage tube, connecting tube, needle, non-return valve.
Capacity: 200ml, 400ml, 600ml, 800ml
Size of Stainless Steel Needle: Fr7, Fr10, Fr12, Fr14, Fr15, Fr16, Fr18, Fr19
Packing: 1pc/Blister bag

200ml	20pcs/43x32x22cm	400ml	12pcs/43x32x22cm
600ml	12pcs/43x32x22cm	800ml	12pcs/54x32x19cm



C033 Urinary Meter Drainage Bag(Double Chambers)

Capacity: 2600ml + 400ml
 With T valve on urine meter
 With twist valve on urine meter
 With needle sample port, with anti-reflux valve
 With reinforced double plastic hanger and rope
 With white bed sheet clamp
Inlet tube: O.D 11 mm, 130cm length
Sterilization: EO
Packing: 1pc/Blister pack, 10pcs/ctn
Carton Size: 48x41x25cm



C031 Urinary Meter Drainage Bag

Capacity:
 Type A 2600ml+500ml, Type B 2000ml+200ml
 With T valve on urine meter, with twist valve on urine meter
 With needle sample port, with anti-reflux valve
 With reinforced double plastic hanger and rope
 With white bed sheet clamp
Inlet tube: O.D 10 mm, 120cm length
Sterilization: EO
Packing: 1pc/Blister pack, 20pcs/ctn
 Type A **Carton Size:** 68x38x28cm
 Type B **Carton Size:** 48x39x31cm



C039 Wound Drainage System

Material: PVC and/or silicone rubber, PP, PS, SS
Application: Used for abdominal, chest, breast and other parts of the fluid, pus and blood drainage.
Components: Container, drainage tube, connecting tube, needle, non-return valve.
Capacity: 200ml, 400ml, 500ml, 800ml.
Size of Stainless Steel Needle: Fr10, Fr12, Fr14, Fr16, Fr18
Packing: 1pc/Blister bag

200ml	20pcs/43x32x22cm	400ml	20pcs/43x32x27cm
500ml	18pcs/48x32x48cm	800ml	12pcs/54x32x19cm



C063 Silicone Wound Drainage System

Material: Silicone rubber, Stainless Steel
Application: Used for abdominal, chest, breast and other parts of the fluid, pus and blood drainage.
Components: Container, drainage tube, connecting tube, needle, non-return valve.
Capacity: 100ml, 150ml, 200ml, 400ml
Drainage Tube: Round type, flat type 7mm & 10mm
Size of Stainless Steel Needle: Fr10- Fr24
Packing: 1pc/Blister bag, 20pcs/box, 80pcs/ctn
Carton Size: 60x41x35cm




C036 Leg Bag

Capacity: 500ml and 1000ml
 Type: ellipse shape
 With green twist valve and green tip cover,
 With anti-reflux valve, with double high elastic belts
Sterilization: EO
Packing: 1pc/PE pack, 200pcs/ctn
Carton Size: 55x40x28cm



C035 Leg Bag

Capacity: 600ml
 Type: Square
 With blue pull-push valve and blue tip cover
 With anti-reflux valve, with double high elastic belts
Sterilization: EO
Packing: 1pc/PE pack, 48pcs/ctn
Carton Size: 30x26x15cm



C064 Chest Drainage Tube

Material: 100% silicone
Feature: Radio-opaque line through the length for X-ray visualization.
 With two holes (round or oval shape) at the end
Size: Fr12- Fr36, length 48.5cm, can be customized.
Packing: 1pc/Blister bag

Fr12-Fr16	40pcs/box	640pcs/ctn	Fr18-Fr24	30pcs/box	480pcs/ctn
Fr26-Fr30	20pcs/box	320pcs/ctn	Fr32-Fr36	15pcs/box	240pcs/ctn

Carton Size: 58x38.5x32.5cm



C074 Drain Thoracic

Material: PVC/TPU
Type: With Trocar, Without Trocar
Pediatric Size: Fr7, Fr8, Fr10, Fr12, Fr14
Adult Size: Fr12- Fr36
Packing: 1pc/blister bag, 10pcs/box
Carton Size: 46.5x34.5x28cm

Ostomy Pouching System

An ostomy pouching system is a prosthetic medical device that provides a means for the collection of waste from a surgically diverted biological system (colon, ileum, bladder) and the creation of a stoma.

Pouch Type: One -piece / Two-piece, Closed or with opening

Base Plate Type: PE foam or Hydrocolloid, with non woven edge or not.

1. Peculiar setuped Hydrocolloid plate has good ability of initial adhesiveness and persistent adhesiveness.
2. It is friendly to skin, no allergic and no irritation. Good resist for PH value between 3 to 12, and resist for different digestive enzyme.
3. Under correct usage, it is safe for long time wearing (5-7 days).
4. Active carbon filter on the bag has the odor removal function.
5. With irrigation port, easy for using and keep the bag clean. With Highly resistant CO-EX film.
6. Flange two piece clasp, can be embedded firmly and tightly, also easy to be separated. Soft, comfortable and easy to use.



C040



C041



C042



C044



C045



C046

C040 One-Piece Closed Pouch

PE foam adhesive base plate
Designed for one time use. Simply discard once full
Dia: 50mm, 60mm
Packing: 20pcs/Box, 20boxes/ctn
Carton Size: 33x31x43cm

C041 One-Piece Closed Pouch

Hydrocolloid adhesive base plate
Designed for one time use. Simply discard once full
Dia: 50mm, 60mm
Packing: 20pcs/Box, 20boxes/ctn
Carton Size: 33x31x43cm

C044 One-Piece Drainable Pouch With Clamp Closure

Hydrocolloid adhesive base plate
With active carbon filter
Type: With irrigation port or without,
With non woven edge or without
Dia: 50mm, 60mm, 70mm, 80mm
Packing: 20pcs/Box, 20boxes/ctn
Carton Size: 33x31x43cm

C042 One-Piece Mini Closed Pouch


Hydrocolloid adhesive base plate
Dia: 50mm, 60mm
Packing: 20pcs/Box, 20boxes/ctn
Carton Size: 33x31x43cm

C045 One-Piece Drainable Pouch With Velcro Closure

Hydrocolloid adhesive base plate
With active carbon filter
Type: With irrigation port or without, with non
woven edge or without
Dia: 50mm, 60mm, 70mm, 80mm
Packing: 20pcs/Box, 20boxes/ctn
Carton Size: 33x31x43cm

C046 One-Piece Drainable Pouch With Twist-Tie Closure

Hydrocolloid adhesive base plate
With active carbon filter
Type: With irrigation port or without, with non
woven edge or without
Dia: 50mm, 60mm, 70mm, 80mm
Packing: 20pcs/Box, 20boxes/ctn
Carton Size: 33x31x43cm



C049 One-Piece Drainable Pouch (Child)

Hydrocolloid adhesive base plate
It is for pediatric stoma.
Dia: 50mm
Packing: 20pcs/Box, 20boxes/ctn
Carton Size: 33x31x43cm



C050 Two-Piece Closure Pouch

Embedded, bag only
Soft and comfortable, Embedded ring. In
miniature size and easy to use.
Dia: 50mm, 60mm
Packing: 20pcs/Box, 20boxes/ctn
Carton Size: 33x31x43cm



C052 Two-Piece Drainable Pouch With Velcro Closure

Embedded, bag only
With active carbon filter
Dia: 50mm, 60mm, 70mm
Packing: 20pcs/Box, 20boxes/ctn
Carton Size: 33x31x43cm



C053 Two-Piece Drainable Pouch With Clamp Closure

Embedded, bag only
With active carbon filter
Dia: 50mm, 60mm, 70mm
Packing: 20pcs/Box, 20boxes/ctn
Carton Size: 33x31x43cm



C054 Two-Piece Drainable Pouch With Twist-Tie

Embedded, bag only
With active carbon filter
Dia: 50mm, 60mm, 70mm
Packing: 20pcs/Box, 20boxes/ctn
Carton Size: 33x31x43cm



C055 Two-Piece System Hydrocolloid Base Plate

With non woven edge
Dia: 50mm, 60mm, 70mm
Packing: 5pcs/box, 20boxes/ctn
Carton Size: 39x23x27cm



C056 Two-Piece System Hydrocolloid Base Plate

Without Non woven edge
Dia: 50mm, 60mm, 70mm
Packing: 5pcs/Box, 20boxes/ctn
Carton Size: 39x23x27cm